Gene Company Limited
基因有限公司 A Gene Group Company

Lonza Nucleofector™ 4D核转染系统

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- **Experiments** – display user-defined experiments
- **Custom programs** – create and manage custom Nucleofection™ Programs
- **Results** – reopen or load result files onto a USB stick plugged into the USB port in front of the Core Unit
- **Lonza programs** – version information about the cell type list
- **Display and audio** – adjust brightness of the touch screen display, time to activate standby and deactivate the touch tone
- **Date and time** – set and format date and time
- **Language and keyboard** – change language settings for display functions and keyboard
- ... more: To switch between pages of the settings list
- **Firmware update** – update the system software (see chapter 2.10 for details)
- **Device cleaning** – opens the tray of the functional modules to remove it for cleaning
- **System restart** – restart system
- **Transport lock** – activates the transportation lock of the system
- **Synchronize** – supports synchronization of the actual firmware with the PC Editor
- **Version** – displays firmware version and serial number of all Modules

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问题	可能原因	建议
Low survival rate	Cells were kept in Nucleofector™ Solution too long	Transfer cells immediately into pre-warmed medium as recommended in the optimized protocol.
	Cells were damaged by harvesting procedure or through handling	Avoid harsh conditions during cell harvesting, especially centrifugation at higher speed or

		overexposure to trypsin. Pipette cells smoothly as they are quite stressed already. Use a plastic Pipette as recommended in the optimized protocols.
	Cells culture conditions were suboptimal	Cells should be viable and in culture for several passages. Avoid excessive cell densities or cell confluencies since this may decrease cell viability post Nucleofection™. For further details please refer to the dedicated optimized protocol.
	Multiple use of Cuvettes	We strongly recommend using the Nucleofection™ Vessels only once, because the high voltage pulses that are applied drastically affect their physical integrity.
	Poor DNA quality	DNA used for Nucleofection™ should be of high purity. We strongly recommend endotoxin-free preparation of the DNA. Do not use procedures involving phenol /chloroform treatment.
Low substrate transfer	DNA amount is too low	We recommend a certain DNA amount per sample (depending on cell type and Nucleofection™ Vessel; for details please refer to respective optimized protocol). If both transfer efficiency and cell mortality are low, the DNA amount could be increased. Increasing DNA amount may lead to higher transfection efficiency, but at the same time result in higher cell mortality.
	Cell number in Nucleofection™	Please use the cell numbers

	sample too high or too low	recommended in the dedicated optimized protocol.
	Poor DNA quality	DNA used for Nucleofection™ should be of high purity. We strongly recommend endotoxin-free preparation of the DNA. Do not use procedures involving phenol /chloroform treatment.

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	Parameter file	Result file, OK	Result file, warning	Result file, error	Skipped well
Undefined					
Defined					
No program control					
No substrate control					