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
The **Object Results** section allows you to save also **results** per **object** to the database. This can be configured individually for each population (local or global). After the evaluation they are available for inspection on the **Evaluation** screen and/or can be exported for use in statistical programs via **Settings – Export Data**. See also section [“Export Data”](#).

#### First level parameter

##### Population

- **None:** **Object results** are not saved in the database to save disk space (default setting).
- **ALL:** All **object** properties which are explicitly calculated by building blocks (e.g. *Calculate Intensity Properties*, *Calculate Morphology Properties*) are saved in the database. I.e. all columns of the **object** tables displayed by the **Define Results** building block are saved.

##### Notice for “ALL”

This option generates large and space-consuming tables in the database. Only choose this option if you really need it.

- **Use Selected Well Results:** Only those cell properties which have corresponding well **results** defined are saved in the database. E.g. if the “mean nuclear area” is defined as a well result, the “nuclear area” is also saved for each individual cell.

##### Notice for “Use Selected Well Results”:

- Properties referenced in formula outputs are not saved in the database.
- The **object** table displayed by the **Define Results** building block shows all properties (not only the selected ones).

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